



Published on *cisTEM* (<https://cistem.org>)

[Home](#) > export Assets and import in Relion

export Assets and import in Relion

Mon, 02/04/2019 - 08:30

#1

rddey80

export Assets and import in Relion

Dear Members,

I exported an asset in relion format from cisTEM and it created a stack file and a star file. Can someone tell step by step how I can use these two files to run 2D classification in Relion? I tried, but it did not work for me. What I did in Relion, is as follows:

1. Imported motion corrected mrc files, that created micrographs.star file
2. used particle extraction as attached screenshot

where to use the stack file?

Thanks for your time,

raj

File:

 [Screenshot1.png](#)

timgrant

Hi Raj,

Hi Raj,

You should be able to straight to the 2D classification job line in Relion, select the star file that you exported from cisTEM, setup your parameters and run.

Cheers!

Tim

Mon, 02/04/2019 - 12:55 [\(Reply to #2\)](#)

rddey80

export Assets and import in Relion

Hi Tim,

Man thanks for your reply. I tried and got an error as follows:

Error:

cannot read file*-stack it does not exist

File: /usr/local/src/relion/src/image.h line: 157

Any idea?

Thanks,

Raj

timgrant

Did you move the stack after

Did you move the stack after you exported it?

Cheers,

Tim

rddey80

Hi Tim,

Hi Tim,

No, I did not move. Both star and stack files are located at the same place where cisTEM exported. Please see below:

```
[rdey@phoebe cistem]$ pwd  
/ddn/TarkD/prex2/cistem  
[rdey@phoebe cistem]$ ls
```

```
micrographs-neg micrographs-pp prex2-225k-stack prex2-225k-  
star.star prex2-neg prex2-pp  
[rdey@phoebe cistem]$ cd ../relion/
```

```
[rdey@phoebe relion]$
```

I ran Relion 2D Classification here browsing the *.star file in the input gui and got the following error

```
ERROR:  
Cannot read file prex2-225k-stack It does not exist  
File: /usr/local/src/relion/src/image.h line: 157
```

Thanks,

Raj

timgrant

Hi Rak,

Hi Raj,

It is strange it does not have an extension, it should have the .mrccs extension. Can you try to export it again with the filename prex2-225k-stack.mrccs

Thanks,

Tim

Tue, 02/05/2019 - 09:59 (Reply to #6)

rddey80

Hi Tim,

Hi Tim,

Thanks. It is now running.

Best regards,

Raj

rddey80

using earlier 2D class average to re-pick particles

Hi Tim,

I have a good ab-initio 3D map generated from a refinement package, which is created from a number of selected 2D class averages (all done in cisTEM). Is it possible to use the 3D ab-initio map or the 2D class averages to re-pick particles in cisTEM? If not, can I switch to relion to do the same? In that case, could you please guide how to do that? I mean, which file to export?

Thanks for your time and help,

Raj

timgrant

Hi Raj,

Hi Raj,

cisTEM does not use references to do the picking, so you can't do this.

Is there a reason why you are keen to do this? Did the picking not work well?

Tim

[Log in](#) or [register](#) to post comments

Source URL: <https://cistem.org/export-assets-and-import-reliion?page=0>