



Published on *cisTEM* (<https://cistem.org>)

[Home](#) > Exporting the list of Good Images

Exporting the list of Good Images

Sat, 01/25/2020 - 12:32

#1

Stephen

Exporting the list of Good Images

Hi Tim,

After Find CTF, I examined the images and the CTF Fit, and selected a set of 'good' images.

Is there a way to export the list of these good images that I have selected.

I would like to also try processing in Relion, but do not want to go through all the images again to select the good ones.

Thanks,

Stephen

timgrant

Hi Stephen,

Hi Stephen,

You can do this by directly accessing the database. This is one of the examples covered in the "images in a group" section here :-

<https://cistem.org/documentation#tab-1-16>

Thanks,

Tim

Sat, 01/25/2020 - 16:19 [Reply to #2](#)

Stephen

Awesome! Thanks!!!

Awesome! Thanks!!!

Tue, 03/10/2020 - 09:30 (Reply to #3)

priyanka

Hi Tim,

Hi Tim,

I am not able to export the good images selected after CTF fit curve. Can you please help for the same.

timgrant

Hi Priyanka,

Hi Priyanka,

Did you try the database commands in the link above? What problem are you seeing?

Thanks,

Tim

Thu, 04/02/2020 - 08:18 [\(Reply to #5\)](#)

priyanka

reply#5

Hi Tim

Initially I was not aware of how to use it. But now problem solved.

Thankyou for your help and concern

Tim

priyanka

regarding import of relion stack into cisTEM

Hi

Can you please help me how to import the stack from relion to cisTEM. When I am trying to do from refinement package import option, I am facing an error.

Thanks

priyanka

Thu, 04/02/2020 - 16:55 (Reply to #7)

timgrant

Hi Priyanka,

Hi Priyanka,

Is this from Relion 3.1? If so then it won't work as 3.1 is not supported in the released version of cisTEM.

There is some discussion here :-

<https://cistem.org/cant-import-relion-31-particles-cistem>

You can modify the star file so that it will import, you essentially have to convert it into the old style.

Thanks,

Tim

Fri, 04/03/2020 - 04:09 (Reply to #8)

priyanka

Hi Tim

Hi Tim

Yes, it is from Relion 3. I will try to modify the file.

Thankyou

Tim

priyanka

regarding alignment of movie file

Hi

Can someone please suggest me regarding an issue with movie alignment. I have encountered a difference between the motion corrected micrograph when I am doing alignment of movies in Relion and cisTEM. I have seen that when I am doing motion correction as well as ctf correction both in cisTEM then resolution fit curve is different as compared to when I am using motion corrected micrographs by Relion then doing ctf correction in cisTEM. I know in both the software the algorithm is different so can someone tell me which way I should follow.

It would be really help in data processing.

Thanks

Priyanka

Tue, 04/14/2020 - 15:04 (Reply to #10)

timgrant

Hi Priyanka,

Hi Priyanka,

The alignments will be different, also in cisTEM the default is to use the movies for CTF estimation, if you import images from Relion then it cannot use Movies.

Thanks,

Tim

priyanka

Reply#8

Hi Tim

I am using Relion 3.0 not Relion 3.1 still I can't import the stack file .Can you please comment on this.

Thanks

Tue, 04/14/2020 - 15:04 (Reply to #12)

timgrant

Hi Priyanka,

Hi Priyanka,

What error are you seeing?

Can you paste the header of the star file here?

Tim

priyanka

reply#13

Hi Tim

This is the error

```
home/grigoriefflab/compiled/wxWidgets3_static/include/wx-  
3.0/wx/strvararg.h(456): assert "(argtype &  
(wxFormatStringSpecifier<T>::value)) == argtype" failed in  
wxArgNormalizer(): format specifier doesn't match argument type.
```

Please suggest me how to tackle it

priyanka

reply#11

Hi Tim

the alignment would be different but in CTF estimation it can take images also. So we can use Relion motioncorrected images.

Wed, 04/15/2020 - 11:58 [\(Reply to #15\)](#)

timgrant

Yes, this will work - I'm

Yes, this will work - I'm just saying why the results may be different.

Tim

timgrant

Can you post the header of

Can you post the header of your star file?

Tim

priyanka

#17

Hi Tim

Sorry to ask but can you please tell what exactly you want me to post

Priyanka

Wed, 04/15/2020 - 13:55 [\(Reply to #18\)](#)

timgrant

The first 100 lines or so of

The first 100 lines or so of your star file that you are trying to import..

Thanks,

Tim

priyanka

regarding retracing the particles selected from 2d classes

Hi

Can someone suggest me how to retrace the particle position selected after 2d classification. I want to see the selected particles on micrograph itself from which particle belongs to.

Thankyou

Priyanka

Sun, 04/19/2020 - 03:23 [\(Reply to #20\)](#)

priyanka

reply#20

Can someone reply on this.It would be helpful

priyanka

reply#19

Hi Tim

This is star file header.

RELION; version 3.0

data_

loop_

_rlnCoordinateX #1

_rlnCoordinateY #2

_rlnClassNumber #3

_rlnAutopickFigureOfMerit #4

_rlnAnglePsi #5

_rlnImageName #6

_rlnMicrographName #7

_rlnMagnification #8

_rlnDetectorPixelSize #9

_rlnCtfMaxResolution #10

_rlnCtfFigureOfMerit #11

_rlnVoltage #12

_rlnDefocusU #13

_rlnDefocusV #14

_rlnDefocusAngle #15

_rlnSphericalAberration #16

_rlnCtfBfactor #17

_rlnCtfScalefactor #18

_rlnPhaseShift #19

_rlnAmplitudeContrast #20

_rlnOriginX #21

_rlnOriginY #22

_rlnGroupNumber #23

_rlnAngleRot #24

_rlnAngleTilt #25

_rlnNormCorrection #26

_rlnLogLikeliContribution #27

_rlnMaxValueProbDistribution #28

_rlnNrOfSignificantSamples #29

_rlnGroupName #30

1963.926503 628.969926 1 1.748691 87.786941 000049@Extract

/job027/Stack/20200304_1125_A013_G000_H006_D001.mrcs

MotionCorr/job010/Stack/20200304_1125_A013_G000_H006_D001.mrc

10000.000000 1.200000 6.523637 -0.00381 200.000000 24168.191406

22763.585938 11.292921 2.700000 0.000000 1.000000 0.000000

0.100000 -2.72922 -10.72922 1 42.809557 107.624284 0.408273

8612.520438 0.751557 3 group_1

Priyanka

Fri, 04/17/2020 - 15:26 (Reply to #22)

timgrant

Hi Priyanka,

Hi Priyanka,

This looks like a star file for images, not a particle stack. Am I correct. cisTEM will only directly import a particle stack and corresponding star file.

Tim

priyanka

reply#22

Hi Tim

This is the particles.star file of selected classes after 2d classification.

Mon, 04/20/2020 - 12:47 [\(Reply to #24\)](#)

timgrant

Hi Priyanka,

Hi Priyanka,

Do you have all your particles in a single file, or do you have a different file for the particles from each image?

Tim

Tue, 04/21/2020 - 05:01 (Reply to #25)

priyanka

Reply#25

Hi Tim

All particles are in one file only.

Priyanka

[Log in](#) or [register](#) to post comments

Source URL: <https://cistem.org/exporting-list-good-images?page=0>