

Home > Exporting or accessing Particle Positions

Exporting or accessing Particle Positions

Wed, 03/14/2018 - 17:03 #1

rasmuskf

Exporting or accessing Particle Positions

Dear cisTEM team,

Good job on this really nice software! I was wondering if there is an option for exporting particle coordinates or somehow accessing the particle positions listed in the table in Assets?

Best,

Rasmus

Re: Exporting or accessing Particle Positions

Hi Rasmus

I've tried the "export" function in cisTEM, writing out particle coords in relion format. It's working well.

"Assets" --> "Refine Pkgs" --> "Export", you could export to relion or frealign format.

 HTH

Ming

Thu, 03/15/2018 - 06:08 (Reply to #2)

rasmuskf

| Hi Ming, |
|---|
| Hi Ming, |
| Thanks for the tip! Don't know why I didn't think of that But thanks again. |
| Best, |

Rasmus

Log in or register to post comments

Source URL: https://cistem.org/exporting-or-accessing-particle-positions