



Published on *cisTEM* (<https://cistem.org>)

[Home](#) > Error: mode 32645 MRC files not currently supported

Error: mode 32645 MRC files not currently supported

Wed, 04/11/2018 - 05:18

#1

bonniemurphy

Error: mode 32645 MRC files not currently supported

Hi,

I'm (recently) having difficulty running CTF estimation on movies. The movies work fine for motion correction, and CTF estimation works fine on the sums that come out of the motion correction, but if I try to run CTF estimation on aligned frames, I get errors like this:

Error: mode 32592 MRC files not currently supported
Error: mode 32747 MRC files not currently supported
Error: mode 32617 MRC files not currently supported
Error: mode 32566 MRC files not currently supported
Error: mode 32567 MRC files not currently supported
Error: mode 32682 MRC files not currently supported
Error: mode 32646 MRC files not currently supported
Error: mode 32580 MRC files not currently supported
Error: mode 32550 MRC files not currently supported
Error: mode 32624 MRC files not currently supported
Error: mode 32528 MRC files not currently supported
Error: mode 32548 MRC files not currently supported
Error: mode 32603 MRC files not currently supported
Error: mode 32688 MRC files not currently supported
Error: mode 32530 MRC files not currently supported
Error: mode 32608 MRC files not currently supported
Error: mode 32713 MRC files not currently supported
Error: mode 32762 MRC files not currently supported
Error: mode 32538 MRC files not currently supported
Error: mode 32712 MRC files not currently supported
Error: mode 32526 MRC files not currently supported
Error: mode 32656 MRC files not currently supported

I have previously been able to run CTF estimation on similar micrographs (ie. same microscope/detector/mag), but I have recently added new movies to the dataset and it's now throwing errors. As an aside, I'm also not able to select a particular micrograph group and run ctf estimation on only the new micrographs - although I've created groups of micrographs that I can select elsewhere in the gui (eg. under actions -> align movies I can select the group 20180404, but under actions -> find ctf I only have the option to select All Images).

Any thoughts on what's going wrong here?

Thank you,

Bonnie

timgrant

Hi Bonnie,

Hi Bonnie,

That error probably means that the program cannot find the movies, is it possible they have been deleted or moved? Is this running locally or on a cluster? If it is on a cluster, are the movies accessible by all nodes?

There is a bug in the current version, that means sometimes you cannot run movie analysis on sub groups - this should be fixed in the next version.

Thanks,

Tim

Mon, 04/16/2018 - 03:42 [\(Reply to #2\)](#)

bonniemurphy

Hi Tim,

Hi Tim,

Thanks for your reply - as you said, it was a problem accessing the files.

Cheers,

Bonnie

[Log in](#) or [register](#) to post comments

Source URL: <https://cistem.org/error-mode-32645-mrc-files-not-currently-supported>